

Go!Link

Fact Sheet

Description: Go!Link is a 12 bit, single-channel interface that connects directly to a USB port on a Windows or Macintosh computer.

Key Features: Compatible with 36 Vernier sensors:

Temperature	Exercise Heart Rate	Flow Rate
Light	Gas Pressure	Force Plate
pH	Magnetic Field	Instrumentation Amplifier
Barometer	Oxygen Gas	Ion-Selective Electrodes
Accelerometer	Relative Humidity	ORP Sensor
Conductivity	Sound Level	Surface Temperature
Current	Turbidity	Extra-Long Temperature
Voltage	UVA Light	Thermocouple
Dissolved Oxygen	UVB Light	Respiration Monitor Belt
Force	Colorimeter	
EKG	Electrode Amplifier	

Specifications Maximum number of sample readings per second: 200
May also be used with Logger Pro 3 software

Audience: Ages 7+

Price: \$59
Includes Logger Lite software (includes a school site license and a student site license for home use)

Installation: Go!Link and Logger Lite come with a Quick Start Installation Guide

Company: Vernier Software & Technology